Pelvic Clock® 4-WEEK EXERCISE PROGRAM



The Pelvic Clock® Exercise Device is for people dedicated to pursuing a regular exercise program to treat their chronic lower back pain and other related conditions.

Some people feel the effects of using the Pelvic Clock® device right away, especially after a long day of sitting. However, please keep in mind that it is not a one-off effort and in order to get the full benefits, you need to make pelvic clock exercises part of your daily regimen.

To help you with that goal, we're providing a basic Pelvic Clock® Exercise Program to allow you to build up your skills using the device and gradually ease into a regular exercise schedule.

You'll find that sticking to this 4-week regimen will not only help you pace yourself, but it will also allow you to see how your body is responding to these pain-relieving stretches over time.

TIPS:

- 1. Start at a slow pace. Take time to build up your exercise routine.
- 2. Listen to your body. When you are looking for pain relief, your own intuition and common sense can be very helpful. You need to find the direction of the movement that reduces your pain.

If tilting your pelvis in one direction increases your pain, tilt it in the opposite direction.

- 3. If you are experiencing pain, stop your practice. Ice your lower back for 15 minutes, and give it a little rest.
- 4. Always exercise within your own limits. Gradually increase your range of motion until you are pain free.
- 5. In the first 2 weeks of your practice, you might want to use a small towel on top of the device for extra padding.
- 6. If your mobility is limited, start exercising on the bed or massage table. The soft surface will reduce the height of the device.


Exercises are a long-term solution.

Work hard and be patient.


How to Set Up For Exercises


Lie on your back with your knees bent and feet flat on the mat (this is called the semi-supine position). The feet should be hip distance apart.

Lift your pelvis about 1 to 3 inches off the floor. Then, place the device below your sacrum, between your sacrum and the floor. Lower your pelvis, and your sacrum should be sitting in the recess of the flat upper surface. The convex lower surface of the device stays in contact with the floor. Point the 12 o'clock marker on the device toward your head.




Once you have the device in position, you can just relax on it for a minute or two. Stretch out your legs on the mat or try rocking side to side.





Exercise Patterns

In each exercise, your pelvis will follow a different pattern, using the numbers of the clock as a guide.




Duration of Exercises

Aim to repeat each exercise 10 to 20 times. The number of repetitions and length of practice will vary for each person, depending on factors like affected side, physical condition, age, gender, and degree of pain. The goal is to perform all of these exercises within a range of motion that does not increase your pain and with the expectation that the device will help you to gradually increase your range of motion until you are pain free.

Find Neutral Pelvis Position

Before you start with Pelvic Clock® Core Stability exercises, you should find your Neutral Pelvis Position when your pelvis is neither rotated nor tilted forward, back, or sideways.

Take five deep breaths, as you balance in this position, keeping the flat surface of the device perfectly parallel to the floor.



Week 1

Knees-to-Chest Stretch Knee-to-Chest Stretch Exercise "3-9" with straight legs Exercise "12-6"

Knees-to-Chest Stretch

Pull your knees toward your chest.

Slowly rock your pelvis in different directions for 30+ seconds.




Knee-to-Chest Stretch

Pull your **right** knee toward your chest. Stretch your left leg and relax.


Slowly rock your pelvis in different directions for 30+ seconds.

Then repeat this stretch with your left knee.



Exercise "3-9" with straight legs

Stretch your legs, relax, and rock your pelvis side to side following the **3-9** pattern.




Exercise "12-6"

Bend your knees and press your feet into the ground.

Slowly tilt your pelvis back toward the **12** o'clock marker (in the direction of your head).

Then take a deep breath and tilt your pelvis forward toward the 6 o'clock marker (in the direction of your feet).



Week 2

Knees-to-Chest Stretch
Knee-to-Chest Stretch
Exercise "3-9" with straight legs
Exercise "12-6"
Exercise "6-1-6-11"
Exercise "Arms Up"
Exercise "Balance on Heels"

Knees-to-Chest Stretch

Pull your knees toward your chest.

Slowly rock your pelvis in different directions for 30+ seconds.




Knee-to-Chest Stretch

Pull your **right** knee toward your chest. Stretch your left leg and relax.


Slowly rock your pelvis in different directions for 30+ seconds.

Then repeat this stretch with your left knee.



Exercise "3-9" with straight legs

Stretch your legs, relax, and rock your pelvis side to side following the **3-9** pattern.




Exercise "12-6"

Bend your knees and press your feet into the ground.

Slowly tilt your pelvis back toward the **12** o'clock marker (in the direction of your head).

Then take a deep breath and tilt your pelvis forward toward the 6 o'clock marker (in the direction of your feet).




Exercise "6-1-6-11"

Slowly rock forward and back following 6-1-6-11 pattern.

6 -1 Tilt forward toward the **6** o'clock marker, then lean back toward **1** o'clock (a back-left tilt).

6- 11 Tilt forward toward the **6** o'clock marker, then lean back toward **11** o'clock (a back-right tilt).



Arms Up

Stretch your arms up. Hands aligned with shoulders.

Take five deep breaths, as you balance in this position, <u>keeping</u> the flat surface of the device perfectly parallel to the floor.



Balance on Heels

Stretch your arms up. Hands aligned with shoulders.

Lift the soles of your feet up and try to balance on your heels.

Take five deep breaths, as you balance in this position, keeping the flat surface of the device perfectly parallel to the floor.



Week 3

Knees-to-Chest Stretch
Knee-to-Chest Stretch
Exercise "3-9" with straight legs
Exercise "12-6"
Exercise "6-1-6-11"
Exercise "12-3-6-9" Clockwise
Exercise "12-9-6-3" Counterclockwise
Exercise "Balance on Heels"
Exercise "Balance on One Heel"

Knees-to-Chest Stretch

Pull your knees toward your chest.

Slowly rock your pelvis in different directions for 30+ seconds.




Knee-to-Chest Stretch

Pull your **right** knee toward your chest. Stretch your left leg and relax.


Slowly rock your pelvis in different directions for 30+ seconds.

Then repeat this stretch with your left knee.



Exercise "3-9" with straight legs

Stretch your legs, relax, and rock your pelvis side to side following the **3-9** pattern.




Exercise "12-6"

Bend your knees and press your feet into the ground.

Slowly tilt your pelvis back toward the **12** o'clock marker (in the direction of your head).

Then take a deep breath and tilt your pelvis forward toward the 6 o'clock marker (in the direction of your feet).




Exercise "6-1-6-11"

Slowly rock forward and back following 6-1-6-11 pattern.

6 -1 Tilt forward toward the **6** o'clock marker, then lean back toward **1** o'clock (a back-left tilt).

6- 11 Tilt forward toward the **6** o'clock marker, then lean back toward **11** o'clock (a back-right tilt).



"Clockwise & Counterclockwise"

Clockwise rotation follows the **12-3-6-9** pattern. Counterclockwise rotation follows the **12-9-6-3** pattern.

Imagine that the axis of the rotation extends from your navel through the center of the device. Then move your pelvis and hips in a circular direction around that center of rotation.



Balance on Heels

Stretch your arms up. Hands aligned with shoulders.

Lift the soles of your feet up and try to balance on your heels.

Take five deep breaths, as you balance in this position, keeping the flat surface of the device perfectly parallel to the floor.




Balance on One Heel

Stretch your arms up. Hands aligned with shoulders.

Lift the soles of your feet up and balance on your heels, keeping the flat surface of the device perfectly parallel to the floor.

Lift your right leg up and bend the right knee at a 90 degree angle. Take five deep breaths, as you balance in this position.

Lift your left leg up and bend the left knee at a 90 degree angle. Take five deep breaths, as you balance in this position.




Week 4

Knees-to-Chest Stretch
Knee-to-Chest Stretch
Exercise "3-9" with straight legs
Exercise "12-6"
Exercise "6-1-6-11"
Exercise "12-3-6-9" Clockwise
Exercise "12-9-6-3" Counterclockwise
Exercise "Balance on Heels"
Exercise "Balance on One Heel"
Exercise "Dead Bug"

Knees-to-Chest Stretch

Pull your knees toward your chest.

Slowly rock your pelvis in different directions for 30+ seconds.




Knee-to-Chest Stretch

Pull your **right** knee toward your chest. Stretch your left leg and relax.


Slowly rock your pelvis in different directions for 30+ seconds.

Then repeat this stretch with your left knee.



Exercise "3-9" with straight legs

Stretch your legs, relax, and rock your pelvis side to side following the **3-9** pattern.




Exercise "12-6"

Bend your knees and press your feet into the ground.

Slowly tilt your pelvis back toward the **12** o'clock marker (in the direction of your head).

Then take a deep breath and tilt your pelvis forward toward the 6 o'clock marker (in the direction of your feet).




Exercise "6-1-6-11"

Slowly rock forward and back following **6-1-6-11** pattern.

6 -1 Tilt forward toward the **6** o'clock marker, then lean back toward **1** o'clock (a back-left tilt).

6- 11 Tilt forward toward the **6** o'clock marker, then lean back toward **11** o'clock (a back-right tilt).



"Clockwise & Counterclockwise"

Clockwise rotation follows the **12-3-6-9** pattern. Counterclockwise rotation follows the **12-9-6-3** pattern.

Imagine that the axis of the rotation extends from your navel through the center of the device. Then move your pelvis and hips in a circular direction around that center of rotation.



Balance on Heels

Stretch your arms up. Hands aligned with shoulders.

Lift the soles of your feet up and try to balance on your heels.

Take five deep breaths, as you balance in this position, keeping the flat surface of the device perfectly parallel to the floor.




Balance on One Heel

Stretch your arms up. Hands aligned with shoulders.

Lift the soles of your feet up and balance on your heels, keeping the flat surface of the device perfectly parallel to the floor.

Lift your right leg up and bend the right knee at a 90 degree angle. Take five deep breaths, as you balance in this position.

Lift your left leg up and bend the left knee at a 90 degree angle. Take five deep breaths, as you balance in this position.



Dead Bug

Lift your legs up in the air with your knees bent at a 90 degree angle.

Stretch your arms upwards. Hands aligned with shoulders.

Take five deep breaths, as you balance in this position, keeping the flat surface of the device perfectly parallel to the floor.

